



## Flask, Volumetric, Polymethylpentene, 100 mL

<b>Product #:</b>	<b>2099542</b>
ZAR Price:	Contact Hach
Ships within 3-5 weeks	

Durable polymethylpentene Class A flask with 14/23 tapered Polypropylene stopper. The stopper provides a firm seal to avoid the possibility of external contamination. Flasks may be autoclaved to 121 °C without affecting the graduation tolerance.

---

### Specifications

Body material:	Polymethylpentene
Capacity:	100 mL
Certifications:	Not Applicable
Type:	Volumetric
Units:	1 each