



## Cylinder, graduated, mixing, glass, 25 mL $\pm 0.3$ mL, 0.5 mL divisions, polyethylene stopper

**Product #:** 2088640  
**ZAR Price:** Contact Hach  
**Ships within 2 weeks**

This glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes a polyethylene stopper. Graduations are in durable blue enamel bake.

---

### Specifications

Body material:	Glass
Capacity:	25 mL
Certifications:	Class B
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing