



Cone, Imhoff settling, 1 L, polycarbonate, w/ screw closure

Product #:	206700
ZAR Price:	Contact Hach
Ships within 1 week	

Fully graduated cones for determination of small amounts of settleable matter in sewage by Imhoff Sedimentation Method. Glass cone features reinforced tip for ruggedness, while the Polycarbonate cone features a leak-proof, removable closure for withdrawal of sediment for analysis. Support ring, support stand, and cone holders available and sold separately.

Specifications

Capacity:	1 L
Quantity:	1 each