



Cylinder, graduated, 10 mL +- 0.1 mL, 0.2 mL divisions (white markings)

Product #: 50838
ZAR Price: Contact Hach
Ships within 2 weeks

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Specifications

Body material:	Glass
Capacity:	10 mL
Certifications:	Class B
Height:	135 mm
Measurement method:	To Deliver
Quantity:	1 each
Type:	Single Metric