



## Cylinder, Graduated, Mixing, Glass, 10 mL $\pm 0.1$ mL, 0.2 mL divisions, glass stopper #13

Product #:	2088638
ZAR Price:	Contact Hach
Available	

This glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

---

### Specifications

Body material:	Glass
Capacity:	10 mL
Certifications:	Class B
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing