



Cylinder, graduated, 50 mL +/- 0.4 mL, 1.0 mL divisions (white markings)

Product #: 50841
ZAR Price: Contact Hach
Ships within 1 week

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Specifications

Body material:	Glass
Capacity:	50 mL
Certifications:	Class B
Height:	190 mm
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric