



Be Right™



Flask, Erlenmeyer, Polycarbonate Capacity 250 mL

Product #:	2089746
ZAR Price:	Contact Hach
Ships within 1 week	

Graduated (Molded-in graduations) flask features: clarity, toughness and strength. Autoclavable; however, repeated autoclaving may result in loss of mechanical strength. #6 Hollow stoppers (P/N 14480-06; pk/6) available separately. 250 ml.

Specifications

Capacity:	250 mL
Material:	Polycarbonate
Units:	1 each