



Buret, Automatic, Self-Zeroing, 50 mL Buret Capacity

Product #: 2087641

OBSOLETE ITEM
This item is no longer available.

Durable, crystal-clear acrylic buret is calibrated to meet ASTM E287 requirements. Leak-proof stopcock with self-lubricating PTFE TFE plug. Unaffected by dilute mineral acids and bases (except hydrofluoric acid and ammonium hydroxide). Handles liquids with a pH range of 1-14 and up to 5% wt/wt of acids or 30% bases can be used. Not designed for exposure to alcohol or organic solvents. Includes low density polyethylene squeeze bottle and 36" tubing for easy self-zeroing set up for titration work stations. Dust cover included.

Specifications

Body material:	Glass
Capacity:	50 mL
Special Features:	Buret, Automatic, Self-Zeroing