



Beaker, Polypropylene, Low Form, 600 mL, 4/pk

Product #:	108082
ZAR Price:	Contact Hach
Ships within 1 week	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifications

Body material:	Polypropylene
Capacity:	600 mL
Diameter:	110.7 mm with spout
Material:	Polypropylene, Low Form
Units:	4 /pk