



Cylinder, Graduated, Polypropylene, 250 mL \pm 1.4 mL, 2.0 mL divisions

Product #:	108146
ZAR Price:	Contact Hach
Ships within 1 week	

Polypropylene construction of both cylinder and base. All sizes have generous pour spout, large molded-in, and easy-to-read graduations. Can be chemically sterilized without affecting accuracy. Autoclavable, but accuracy will be adversely affected. Cylinders meet accuracy requirements of ASTM Class B, E1272, ISO Standard 6706, and CFR 21, Part 177.1520

Specifications

Body material:	Polypropylene
Capacity:	250 mL
Certifications:	Class B
I. D.:	34.5 mm
Length:	353 mm
Measurement method:	To Deliver/To Contain
Quantity:	1 each
Type:	Single Metric