



## Cylinder, Graduated, 25 mL, Certified

<b>Product #:</b>	<b>2636240</b>
ZAR Price:	Contact Hach
Ships within 2 weeks	

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes  $\geq$  25 mL.

---

### Specifications

Body material:	Glass
Capacity:	25 mL
Certifications:	Class A
Height:	192 mm
Measurement method:	To Deliver
Quantity:	1 each
Type:	Single Metric