



Bottle, BOD, 300 mL, 24/pk, serialized: 49 - 72

Product #:	2898900
ZAR Price:	Contact Hach
Available	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stoppers included

Serialized bottles have silk-screened, sequential numbers (49-72) and include a bar code.

Specifications

Body Material:	Glass
Capacity:	300 mL
Material:	Borosilicate Glass
Package size UOM:	24/pk, Serialized: 49 - 72
Quantity:	24 /pk