



## Cylinder, graduated, mixing, Glass, 50 mL $\pm$ 0.4 mL, 1.0 mL divisions, glass stopper #16

<b>Product #:</b>	<b>2088641</b>
ZAR Price:	Contact Hach
Ships within 1 week	

This tall form glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

---

### Specifications

Body material:	Glass
Capacity:	50 mL
Certifications:	Class B
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing