



Flask, Volumetric, Class A, Certified, 500 mL

Product #:	2636649
ZAR Price:	Contact Hach
Ships within 2 weeks	

Individually serialized glass flask supplied with Certificate of Graduation Accuracy and traceable to NIST Standards. Class A tolerances meet requirements of ASTM E438 and ASTM E288. Graduation ring on neck of flask. Supplied with ground glass stopper. (replacement stoppers size 19).

Specifications

Body material:	Class A Glass
Capacity:	500 mL
Certifications:	Yes
Units:	1 each