



## Bottle, BOD, 300 mL, 24/pk, serialized: 1-24

<b>Product #:</b>	<b>2898700</b>
ZAR Price:	Contact Hach
Available	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included Serialized bottles have silk-screened, sequential numbers (1-24) and include a bar code.

---

### Specifications

Body Material:	Glass
Capacity:	300 mL
Material:	Borosilicate Glass
Package size UOM:	24/pk, Serialized: 1 - 24
Quantity:	24 /pk