



## Cylinder, graduated, 100 mL +- 0.6 mL, 1.0 mL divisions (white markings)

Product #:	50842
ZAR Price:	Contact Hach
Ships within 1 week	

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes  $\geq$  25 mL.

---

### Specifications

Body material:	Glass
Capacity:	100 mL
Certifications:	Class B
Height:	255 mm
Measurement method:	To Deliver
Quantity:	1 each
Type:	Single Metric