



## Centrifuge Tubes, Round Bottom, 10/pk

Product #:

2278739

 **Hazardous**

**OBSOLETE ITEM**

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Round bottom, polycarbonate centrifuge tube. 15 mL (1.6 cm x 11.4 cm). Working range -135 - 135 °C (-211 - 275 °F). Nontoxic. Autoclavable. Maximum RCF is 50,000 xG. Note: For proper performance, fill tubes to at least 80% of total capacity. No graduations. No Closure included, purchase item number 2585220. Autoclavable.

---

### Specifications

Capacity:	15 mL
Length:	115 mm
Quantity:	10 /pk