



Cylinder, Graduated, Mixing, 10 mL, Certified

Product #:

2636338

 **Hazardous**

OBSOLETE ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Specifications

Body material:	Glass
Capacity:	10 mL
Certifications:	Class A
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing