



## Cylinder, Graduated, Polymethylpentene, 1000 mL $\pm$ 5.0 mL, 10.0 mL divisions

**Product #:** 217253  
**ZAR Price:** Contact Hach  
**Ships within 1 week**

This cylinder features transparent, nontoxic, noncontaminating, and biologically inert material with excellent chemical resistance. Meets ASTM E1272 Class B accuracy requirements, ISO Standard 6706, and suitable for food and beverage use (PMP resin meets requirements of CFR 21, Part 177.1520). Calibrated “To Contain/To Deliver” at 20°C (68°F). Shatterproof and thermally resistant to 135°C (275°F). Absence of meniscus for easy-to-read levels on front-to-back raised graduations.

---

### Specifications

Body material:	Polymethylpentene
Capacity:	1000 mL
Certifications:	Class B
Measurement method:	To Deliver/To Contain
Quantity:	1 each
Type:	Single Metric