



Buchner Funnel 186 mL approx. volume, Polypropylene

Product #: 2090553

OBSOLETE ITEM
This item is no longer available.

Light weight and durable polypropylene construction with a removable filter plate for easy cleaning, chemical resistant. Fully autoclavable. Stopper sold separately.

Capacity: 186 mL

Inside Top Diameter: 76 mm

Outer Diameter of stem: 11 mm

Stem Length: 71 mm

Overall Height: 146 mm

Fits Filter Paper Size: 70 mm

Autoclavable: up to 135 °C

Specifications

Approximate Volume:	186 mL
Body material:	Polypropylene
Capacity:	186 mL
I. D.:	76 mm
Length:	Stem: 65 mm