



## Beaker, Griffin, Low Form, Glass, 1000 mL, 6/pk

<b>Product #:</b>	<b>50083</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Specifications

Body material:	Glass
Capacity:	1000 mL
Material:	Glass, Griffin, Low Form
Units:	6 /pk