



Beaker, Griffin, Low Form, Glass, 400 mL, 12/pk

Product #:	50078
ZAR Price:	Contact Hach
Available	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifications

Body material:	Glass
Capacity:	400 mL
Material:	Glass, Griffin, Low Form
Units:	12 /pk