



Flask, Volumetric, Class A, Certified, 100 mL

| | |
|---------------------|--------------|
| Product #: | 2636642 |
| ZAR Price: | Contact Hach |
| Ships within 1 week | |

Individually serialized glass flask supplied with Certificate of Graduation Accuracy and traceable to NIST Standards. Class A tolerances meet requirements of ASTM E438 and ASTM E288. Graduation ring on neck of flask. Supplied with ground glass stopper. (2603213 size 13 replacement stoppers).

Specifications

| | |
|-----------------|---------------|
| Body material: | Class A Glass |
| Capacity: | 100 mL |
| Certifications: | Yes |
| Units: | 1 each |