



## Cylinder, Graduated, Polycarbonate, 250 mL, 2.0 mL divisions

<b>Product #:</b>	<b>2117946</b>
ZAR Price:	Contact Hach
Available	

Clear polycarbonate cylinder is autoclavable and withstands temperatures to 138°C (280°F). Features easy to read graduations in durable fused ink.

---

### Specifications

Body material:	Polycarbonate
Capacity:	250 mL
Certifications:	None
Measurement method:	To Deliver/To Contain
Quantity:	1 each
Type:	Single Metric