



Cylinder, Graduated, 2000 mL, with Handles

Product #:	2612954
ZAR Price:	Contact Hach
Available	

One-piece polypropylene construction with a molded in handle top ensure a secure, solid grip and superior control for safer, more convenient pouring. The “Total Contained, Total Delivered ”double scale is accurate and easy to read (cylinder will not create a meniscus), while the hexagonal base prevents tipping or rolling. Hand opening is 89x 29 mm (3 x 11/8"). Steam autoclavable upto 121°C (250°F).

Specifications

Body material:	Polypropylene
Capacity:	2000 mL
Certifications:	None
Measurement method:	To Deliver/To Contain
Quantity:	1 each
Type:	Double Metric