



Beaker, Polypropylene, Low Form, 250 mL, 6/pk

Product #:	108076
ZAR Price:	Contact Hach
Available	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifications

Body material:	Polypropylene
Capacity:	250 mL
Diameter:	83.6 mm with spout
Material:	Polypropylene, Low Form
Units:	6 /pk