



## Bottle, BOD, 60 mL capacity, marked 30 mL graduation line

Product #:	190902
ZAR Price:	Contact Hach
Ships within 2 weeks	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 mL sizes).

---

### Specifications

Body Material:	Borosilicate Glass
Capacity:	60 mL
Material:	Borosilicate Glass
Package size UOM:	1pk
Quantity:	1 each