



Cylinder, graduated, mixing, certified, 25 mL

Product #:	2636340
ZAR Price:	Contact Hach
Ships within 1 week	

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Specifications

Body material:	Glass
Capacity:	25 mL
Certifications:	Class A
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing