



Cylinder, Graduated, Polypropylene, 500 mL \pm 2.6 mL, 5.0 mL divisions

| | |
|---------------------|---------------|
| Product #: | 108149 |
| ZAR Price: | Contact Hach |
| Ships within 1 week | |

Polypropylene construction of both cylinder and base. All sizes have generous pour spout, large molded-in, and easy-to-read graduations. Can be chemically sterilized without affecting accuracy. Autoclavable, but accuracy will be adversely affected. Cylinders meet accuracy requirements of ASTM Class B, E1272, ISO Standard 6706, and CFR 21, Part 177.1520

Specifications

| | |
|---------------------|---------------------------|
| Body material: | Polypropylene |
| Capacity: | 500 mL |
| Certifications: | Class B |
| Instrument: | All except DR/4000, PC II |
| Measurement method: | To Deliver/To Contain |
| Method: | 8006 |
| Quantity: | 1 each |
| Type: | Single Metric |