



Flask, Volumetric, Class A, Certified, 250 mL

Product #:	2636646
ZAR Price:	Contact Hach
Ships within 1 week	

Individually serialized glass flask supplied with Certificate of Graduation Accuracy and traceable to NIST Standards. Class A tolerances meet requirements of ASTM E438 and ASTM E288. Graduation ring on neck of flask. Supplied with ground glass stopper. (PN 2603216 size16 replacement stoppers).

Specifications

Body material:	Class A Glass
Capacity:	250 mL
Certifications:	Yes
Units:	1 each