



## Sample Container, 120 mL High PF

<b>Product #:</b>	<b>2552200</b>
ZAR Price:	Contact Hach
Ships within 1 week	

These 120 mL (4oz.) High Profile Corning Snap-Seal™ polypropylene containers are designed for a wide variety of applications and provide a reliable leak-tight seal when closed properly. Approximately 3.5 inches tall (91 mm) and 1.75 inches (45 mm) diameter. Plastic, snap seal, translucent. Sold individually.

---

### Specifications

Capacity:	120 mL - 1 each
Diameter:	45 mm (approx)
Quantity:	1 /pk