



## Cylinder, graduated, mixing, certified, 100 mL

<b>Product #:</b>	<b>2636342</b>
ZAR Price:	Contact Hach
Available	

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

---

### Specifications

Body material:	Glass
Capacity:	100 mL
Certifications:	Class A
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing