



Cylinder, graduated, mixing, certified, 50 mL

Product #:	2636341
ZAR Price:	Contact Hach
Available	

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Specifications

Body material:	Glass
Capacity:	50 mL
Certifications:	Class A
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing