



Disposable BOD bottles, 300 mL, 100/cs

Product #:	2943100
ZAR Price:	Contact Hach
Ships within 1 week	

Disposable BOD Bottle

Case of 100 disposable BOD bottles without stoppers. Plastic bottle is EPA approved for BOD measurement.

Single bottles available under article number 2943101. Plastic stoppers (25 pack) available under article number 2943900. Glass stoppers (quantity of 1) available under article number 190901.

Fresh bottle, ready for use

No cleaning required

Unbreakable and recyclable

Oxygen will not diffuse in or out of sample

EPA compliant

Specifications

Body Material:	Plastic
Capacity:	300 mL
Material:	PET
Method:	8043
Quantity:	100 100/pk