



## Beaker, Polypropylene, Low Form, 1000 mL, 3/pk

<b>Product #:</b>	<b>108083</b>
ZAR Price:	Contact Hach
Available	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes  $\geq 250$  mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Specifications

Body material:	Polypropylene
Capacity:	1000 mL
Diameter:	129.8 mm with spout
Material:	Polypropylene, Low Form
Units:	3 /pk