



Bottle, Sampling, Low Density, Polyethylene, 250 mL, 12/pk

Product #:	2087076
ZAR Price:	Contact Hach
Ships within 1 week	

Low-density polyethylene (LDPE) bottles feature wide-mouth opening for easy filling and leak-proof polypropylene closure. Great choice for grab samples, shipping liquids, soil sampling, and storing dry chemicals.

Specifications

Body Material:	Low Density Plastic
Capacity:	250 mL
Material:	LDPE
Package size UOM:	12pk
Quantity:	12 /pk