



Cylinder, Graduated, 50 mL, Certified

Product #:	2636241
ZAR Price:	Contact Hach
Ships within 1 week	

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in

white marking. Bumper guard included for all sizes ≥ 25 mL.

Specifications

Body material:	Glass
Capacity:	50 mL
Certifications:	Class A
Height:	225 mm
Measurement method:	To Deliver
Quantity:	1 each
Type:	Single Metric