



## Cylinder, Graduated, Polypropylene, 1000 mL $\pm$ 5.0 mL, 10.0 mL divisions

<b>Product #:</b>	<b>108153</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Polypropylene construction of both cylinder and base. All sizes have generous pour spout, large molded-in, and easy-to-read graduations. Can be chemically sterilized without affecting accuracy. Autoclavable, but accuracy will be adversely affected. Cylinders meet accuracy requirements of ASTM Class B, E1272, ISO Standard 6706, and CFR 21, Part 177.1520

---

### Specifications

Body material:	Polypropylene
Capacity:	1000 mL
Certifications:	Class B
Measurement method:	To Deliver/To Contain
Quantity:	1 each
Type:	Single Metric