



Cylinder, graduated, 250 mL +- 1.4 mL, 2.0 mL divisions (white markings)

Product #: 50846
ZAR Price: Contact Hach
Ships within 2 weeks

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Specifications

Body material:	Glass
Capacity:	250 mL
Certifications:	Class B
Height:	330 mm
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric