



## Bottle, BOD, black, 300 mL

Product #:	2086047
ZAR Price:	Contact Hach
Ships within 1 week	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 and 300 mL sizes). Serialized bottles have silk-screened, sequential

---

### Specifications

Body Material:	Glass
Capacity:	300 mL
Material:	Borosilicate Glass
Package size UOM:	1pk
Quantity:	1 each