



3 Place Manifold Kit, Glassware

Product #:

2943232

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Xenosep Manifold Glassware kit for use in conjunction with the 3-Place Stainless Steel Manifold. Manifold not included.

Specifications

What's included?: Kit includes glassware necessary to process three samples at a time when used in conjunction with 2946707 Stainless Steel Manifold.

What's included?

Kit includes glassware necessary to process three samples at a time when used in conjunction with 2946707 Stainless Steel Manifold.