



Replacement: Filtration Assembly for Xenosep SPE Starter Kit

Product #: 2946706

OBSOLETE ITEM
This item is no longer available.

Filtration Assembly includes: 1 x Funnel, 300 mL, 47mm (PN 2943217), 1 x Holder, Filtration, 47mm (PN 2943236), 1 x Stopper, Silicone w/9.5 mm hole, #8 (PN 2943238), 2 x O-Ring, Coupler, 47mm (PN 2943223), 1 x Coupler, PTFE, 47mm (PN 2943218), 1 x Support, Stainless Steel, 47mm (PN 2943219), 1 x Clamp, Aluminum, 47mm, Wide Opening (Pn 2943220).