



Be Right™



pH paper, 9.0 - 12.0 pH range

| | |
|---------------------|--------------|
| Product #: | 38533 |
| ZAR Price: | Contact Hach |
| Ships within 1 week | |

Provides simple, fast, economical method of pH measurement. Sharp color change for each 0.5 pH unit. Convenient for sample pH adjustments. Recommended for sewage water, not finished water. Five 4.6-meter (15') rolls per package.

Specifications

| | |
|-----------|-------------------------------------|
| Quantity: | 5 4.6-meter (15') rolls per package |
| Range: | 9.0 - 12.0 pH units |