



pH Test Strip, 7.5 - 14 pH units, 100 tests

Product #:	2751000
ZAR Price:	Contact Hach
Available	

Most natural waters range from pH 4 - pH 9, but commonly are above pH 7 because of carbonates and bicarbonates (alkalinity). The pH of water used in industry, boiler feed water, and swimming pools usually is kept within a narrow range.

Test Strips are great for quick evaluations of pH levels in the lab or in the field

Easy to use, disposable, and inexpensive!

Specifications

Number of tests:	100
Parameter:	pH
Quantity:	100 /pk
Range:	7.5 - 14 pH units
Smallest Increments Steps:	7.5, 8, 8.5, 9, 9.5, 10, 10.5, 11, 11.5, 12, 12.5, 13, 13.5, 14,
Units:	100 tests