



Xenosep EPA 1664A Standard, 1x30 mL, 4mg/mL

Product #:

2943701

ZAR Price:

Contact Hach

No availability shown

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Xenosep EPA 1664A Standard, 1x30 mL, 4mg/mL. Pipet 10 mL to achieve the 40mg PAR (Precision and Recovery) Standard needed for the method. Each bottle contains enough standard to QC a batch of samples.

Specifications

Quantity: 1 each