



pH Paper, 0 - 14 pH Range, 100 tests

Product #: 2601300
ZAR Price: Contact Hach
Ships within 2 weeks

Provides a simple, fast, economical method of pH measurement. Obtain accurate results in weakly buffered, turbid and lightly colored samples. Measures 0 - 14 pH.

Test Strips are great for quick evaluations of pH levels in the lab or in the field

Easy to use, disposable, and inexpensive!

Specifications

Method:	10056
Number of tests:	100
Parameter:	pH
Platform:	Test Strip
Range:	0 - 14 pH units
Smallest Increments Steps:	0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14
Units:	100 tests