



## Distillation apparatus set

**Product #:** 2265300  
**ZAR Price:** Contact Hach  
**Ships within 2 weeks**

### Fast set-up for sample pretreatment.

Fast set-up for sample pretreatment. This Distillation set includes EPA approved glassware ONLY. Easy to set up because you only receive the parts you need. Illustrated, step-by-step procedures simplify distillation.

\*Heater and support not included.

**Distillation set includes all the glassware you need to satisfy EPA distillation requirements**

**Easy to set up - you only receive the parts you need**

**Illustrated, step-by-step procedures simplify distillation**

---

### Specifications

Quantity:	Each 1
What's included?:	Distillation flask, J-tube, condenser, drip tube, erlenmeyer flask, and non-Hg thermometer

---

### What's included?

Distillation flask, J-tube, condenser, drip tube, erlenmeyer flask, and non-Hg thermometer