



Be Right™



Bottle, Wash, Nalgene, Narrow Mouth, 500 mL, Yellow Cap/ Stem, 6/pk

Product #: 2928201

OBSOLETE ITEM
This item is no longer available.

Narrow mouth polypropylene screw closure and stem feature one-piece construction to ensure a leak proof seal. Tip can be cut back to increase flow. 500 mL (16.8 oz. Capacity, 28 mm mouth). Choose from 4 cap/stem colors.

Specifications

Capacity:	500 mL
Color:	Clear
Material:	LDPE
Quantity:	6 pk
Stem:	Yellow Cap/Stem
Unit:	6 pk