



Bottle, Dropper, Polystop, 60 mL

Product #:

2929902

 **Hazardous**

OBSOLETE ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

The Polystop Dropper Bottle provides a distinct advantage over other glass dropper bottles. The dropper section easily pulls free for use while the polypropylene stopper remains snugly in the neck of the bottle, which helps prevent chipping the glass dropper section.

Specifications

Body Material:	Glass
Capacity:	60 mL
Material:	Glass
Quantity:	1 each