



Be Right™



Bottle, Dropper, Polystop, 100 mL

Product #: 2929903

OBSOLETE ITEM
This item is no longer available.

The Polystop Dropper Bottle provides a distinct advantage over other glass dropper bottles. The dropper section easily pulls free for use while the polypropylene stopper remains snugly in the neck of the bottle, which helps prevent chipping the glass dropper section.

Specifications

Body Material:	Glass
Capacity:	100 mL
Material:	Glass
Package size UOM:	1pk
Quantity:	1 each